

# **A novel Graphical User Interface for Non-Invasive Brain Stimulations**

Gregor Dederichs  
May 2024

## ■ Introduction

- Temporal Interference
- Current Interface Limitations

## ■ Explored Approaches


- Troubleshooting in MATLAB
- NI-MAX software
- Alternative Python packages

## ■ Results – HummelGUI

- Setup
- Functionalities and features
- Code structure & Adding features, Limitations

# Introduction

- Noninvasive Method for Deep-Brain Stimulation
- Temporally Interfering Electric Fields
  - Low frequency envelopes drive neuron spiking
  - Low amplitude envelopes do not recruit neurons
  - Target areas and structures **specifically**





Adapted from Grossmann et al. 2017

# Temporal Interference

## Parameters allow steerability

- Striatum
- Hippocampus



## Implications for Neurological Diseases

- Alzheimer's disease
- Parkinson's disease
- Stroke
- ...

## Controlled through a Graphical User Interface (GUI)

# User Interface Limitations

- **Voltage Discharge**

- Occurs when using a trigger
- Can be **pain**-causing for patients
- **Impractical** for user – Prevents safe termination of stimulations

- **No flexibility**

- Parameters are fixed
- **No parameter update** possible during stimulations

- **Not open source**

- **Decreased accessibility** and sharing

- **Impractical Blinding**

- **Unblinding** risks and **decreased usability**

# Explored Approaches

**DAQ Theory and MATLAB  
Troubleshooting**


**NI-MAX Software**

**Alternative packages  
for Python**

## ■ Voltage Discharge – Probable Source

- Entire signal is stored in DAQ
- Stop → Simultaneous release
  - Without triggers: Internal control of DAQ → Spike Prevented
  - **With triggers:** Error when closing trigger channels

DAQ freezes: Internal control lost → **Voltage Spike**






- Free software by National Instruments
- Device Test Panel
  - Allows testing of simple features
  - Improper trigger setup → **Voltage Spike**
- Virtual Testing
  - Through NI-MAX
  - Virtual devices – Rapidly test certain features without a complete setup
- Device Setup
  - Device **specifications**
  - Device **status** checks
  - Device **configuration** – Renaming devices

- **Package: nidaqmx**
  - API for interacting with the NI-DAQmx driver
  - Open source → Increased **accessibility**
- **Notably limited on current DAQs**
  - Attempts to hard-code functionalities **failed**
    - Example: **No timing** configuration on **USB-6216**
    - Needs explicit specification in Python
    - Idea: Subdivide signal to be able to update it on the fly
    - Problem: **Accumulates delay**

- **Limitations**

- API **created generally** for all NI DAQs
- **Not all DAQs are compatible** with all functionalities



# HummelGUI

For NI USB-6341


# System Setup

## ■ Connections

- **Identical** to current setup
  - BNC 1: "AO0" to stimulator 1
  - BNC 2: "AO1" to stimulator 2
  - BNC 3: "PFI0" to trigger port
  - USB: DAQ to computer running GUI


## ■ Dependencies

- **Different** from current setup
- Python and related **packages**
  - Numpy
  - PyQt6
  - Pyqtgraph
  - Nidaqmx
  - Pandas
  - Matplotlib



## ■ Intended Use

- Troubleshooting
- Testing Parameters
- Checking Functionality
- Titration



## ■ Selecting Stimulation


- Drop down menu
  - iTBS
  - cTBS
  - TBS control
  - TI



The screenshot shows a software interface with a light gray background. At the top, the text "Stimulation Type" is displayed in a bold, black font. To the right of this text is a small, empty square checkbox followed by the letter "B". Below the title is a white rectangular dropdown menu with a thin gray border. Inside the menu, the text "iTBS" is visible, and a small downward-pointing chevron is on the right side. In the bottom right corner of the interface, the word "Setti" is written in a large, white, sans-serif font.

## ■ Adjusting Parameters

- Value Fields
  - Duration
  - Frequency
  - Currents
- Save Option
  - Avoids unnecessary saves during troubleshooting
- Reset Button



**Stimulation Parameters** ☐ Save Parameters

**Duration Parameters**

Total Duration (seconds): 50

Ramp-up Time (seconds): 5

Ramp-down Time (seconds): 5

**Frequency Parameters**

Carrier Frequency (Hz): 2000

Pulse Frequency (Hz): 100

Burst Frequency (Hz): 5

**Current Parameters**

Amplitude Sum (mA): 4

Amplitude Ratio (A1/A2): 1


Amplitude Values: A1 = 2.0 [mA]; A2 = 2.0 [mA]

**Repetition Parameters**

Number of Stimuli:

## ■ Checking Waveform


- Created Through a Button
- **Necessary** before Running Stimulation in **Settings Mode**
- Yellow lines mark signal between ramp-up and –down
- Zoom on x-axis only for easy visualisation





## ■ Running Stimulation

- Choose Number of Repetitions
- Trigger Option
- While Running
  - **Update** stimulation
    - Does not restart repetition count
  - **Stop** stimulation
    - Instantly updates to a ramp down signal
    - **No Discharge**



- **Intended Use**
  - Study Experiments
    - Removes access to parameters and waveforms
    - **No unblinding risk**

**Stimulation Type** ☒ Blind Mode

Select Subject ID  Select Session ID

Get Stimulation from file

**Stimulation Parameters** ☒ Save Parameters

Default Folder to Save Parameters: parameter\_history

Run Stimulation

Update Stimulation Stop Stimulation

☒ Trigger

Check session and subject ID

## ■ Selecting IDs

- Only from Existing Data in excel file
- Load stimulation with button – Check Label

## ■ Saving Options

- Parameters are saved to *parameter\_history* (directory) – Modifiable

The screenshot displays the 'Hummel Lab Graphical User Interface' with the following elements:


- Stimulation Type**: A section header with a checkbox labeled 'Blind Mode' that is checked.
- Select Subject ID**: A dropdown menu.
- Select Session ID**: A dropdown menu.
- Get Stimulation from file**: A button.
- Stimulation Parameters**: A section header with a checkbox labeled 'Save Parameters' that is checked.
- Default Folder to Save Parameters: parameter\_history**: A text input field.

## ■ Running Stimulation

- Identical to Settings Mode
- Trigger option
- Run, update and Stop options

## ■ Label

- In Blind and Settings Modes
- Guides User with all important steps (errors, status, ...)



- Implement **GUI function** and **DAQ communication** (stimulation)

File	Function	Modifications
main.py	Launching GUI	None
GUI.py	GUI Layout	Adding Widgets
	Waveform Creation	Adding Signals – following previous implementations as template – see also utility files
GUI_worker.py	DAQ Communication	Rare – Technical Functionalities (ex: changing trigger settings)
	Running Tasks	None

# Code Structure – Utility Files


- Implement **Underlying Structure and Utility**

File	Function	Modifications
util.py	Storing Defaults	Setting/Adding Defaults
fbase.py	Storing Base Functions	Rare
iTBS.py	Creates Signals of Different Stimulation  One File per Stimulation	Adding new Stimulations – using previous files as template – see also GUI.py
cTBS.py		
TBS_ctrl.py		
Tl.py		

## ■ Available on GitHub

- <https://github.com/gdedericchs/HummelGUI>

- **Quick Guide**
- **Installation** Procedure
- **Setup** Instructions
- Detailed GUI **Functionalities**
- File **Descriptions**
- Future **Modifications** Guide



# Thank You